INNOVATIVE TECHNOLOGY WORLDWIDE

DIAPHRAGM PUMPS WITH KNF STABILIZATION SYSTEM

SECTION 60.14



N940.5APE shown with IP 54 AC motor



N940.5APE-B with multi-voltage power supply input (As an option adjustable with potentiometer or by external signal input).



Concept

The powerful diaphragm vacuum pumps N940 series are especially well suited for all applications requiring excellent suction speed with low absolute pressures. A series of technical innovations, including a diaphragm stabilization system, allows the high suction especially in the low-vacuum range.

The vacuum pump N940.5APE-W is available with optional adjustable flow rate; in this case, either an potentiometer or an external controller with an analog signal input will alter the motor speed. These pump version make it possible to adapt the flow rate to the requirements of a specific applications process.

Features

Transferring and evacuation of air and gases

No contamination of the media due to oil-free operation

KNF stabilization system Optimized suction speed, also for low absolute pressures

High level of gas tightness approx. 6×10^{-3} mbar x l/s (not tested in serial production)

Quiet running

Cool and efficient brushless motor

Multi-voltage power supply input by N940.5APE-W

Can operate in any installed position

Areas of use

The N940 series of diaphragm vacuum pumps offer a high level of performance in a compact unit size. Typical applications are in the fields of analysis, chemistry, medicine and production technologies.

The N940 pumps also support turbomolecular systems as roughing pumps.

There are many applications for the N940 series, please contact KNF for application advice.

PERFORMANCE DATA					
Туре	Delivery (I/min)	Vacuum (mbar absolute)	atm. Press.	Pressure (bar g)	Weight (kg)
N940.5APE	50	1.5		0.5	18.9
N940.5APE-W	50	1.5		0.5	16.8

N940.5APE

PERFORMANCE DATA

Type and Order No.	Delivery (I/min) ¹⁾ at atm. pressure	Max. operating pressure (bar g)	Ultimate vacuum (mbar abs .)
N940.5APE	50	0.5	< 1.5
¹⁾ Litre at STP			

MOTOR DATA

Motor type: AC version			
Protection class		IP 54	
Voltage/Frequencie	es (V/Hz)	115/60	
Power P ₁	(W)	250	
Operating current	(A)	3.4	

Motors with other voltages, frequencies and protection classes on request.

Dimensions mm (All dimensional tolerances conform to DIN ISO 2768-1, Tolerance Class V)





MODEL CODES AND MATERIALS

Type and Order No.	Pump head	Diaphragm	Valves
N940.5APE	Aluminium	EPDM	EPDM

N940.5APE

PERFORMANCE DATA

Type and Order No.	Delivery (I/min) ¹⁾ at atm. pressure	Max. operating pressure (bar g)	Ultimate vacuum (mbar abs.)
N 940.5 APE	50	0.5	< 1.5
	¹⁾ Litre at STP		

MODEL CODES AND MATERIALS

Type and Order No.	Pump head	Diaphragm	Valves
N940.5APE	Aluminium	EPDM	EPDM

The flow rate of the N940.5APE series pumps can be adjusted with an optional potentiometer or by external analog signal. Please contact us for further information.

MOTOR DATA

Motor type: brushless DC motor with AC power supply				
Protection class	IP 20			
Voltage (V)	~90-264/50-60			
Power P1 (W)	180			
Operating current (A)	1.85			

Dimensions mm (All dimensional tolerances conform to DIN ISO 2768-1, Tolerance Class V)





PERFORMANCE DATA

Suction pumping speed



Pump down time for 20 litre receiver



Accessories

Addessories				
Description	Order No.	Details		
Silencer	045993	G 3/8		
Adapter for silencer	014757	G 3/8 to G 1/4		
Small flange, stainless steel	048116	G 1/4, DN 16		
Hose connector, Ms	049880	G 1/4, for tube ID 13		
Hose connector, PP	0045293	G 1/4, for tube ID 10		
Sealing for hose connector, Ms	029112			
Adjustable delivery through speed regulation for N 940.5 APE-W	on request	with potentiometer or analog signal input		

Please visit our website at <u>www.knf.com</u>/usa.htm for further information.

KNF NEUBERGER, INC.

2 Black Forest Road Trenton, NJ 08691-1810 609-890-8600 • Fax 609-890-8323

KNF reserves the right to make changes.

SECTION 60.14 (0706)

INNOVATIVE TECHNOLOGY WORLDWIDE

TIPS ON FUNCTION, INSTALLATION AND SERVICE

THE BASIC FUNCTION OF KNF DIAPHRAGM VACUUM PUMPS AND COMPRESSORS

An elastic diaphragm is moved up and down by an eccentric (see illustration). On the down-stroke it draws the air or gas being handled through the inlet valve. On the up-stroke the diaphragm forces the medium through the exhaust valve and out of the head. The compression chamber is hermetically separated from the drive mechanism by the diaphragm. The pumps transfer, evacuate and compress completely oil-free.

Diaphragm pump

INSTALLATION AND OPERATION

- Range of use: Transferring air and gases at temperatures between + 5 $^\circ\mathrm{C}$ and + 40 $^\circ\mathrm{C}$
- Permissible ambient temperature: between + 10 °C and + 40 °C
- Standard pumps are not suitable for use in areas where there is a risk of explosion. In these cases there are other products in the KNF program please ask us for details
- To prevent the maximum operating pressure being exceeded, restriction or regulation of the air flow should only be carried out in the suction line
- Components connected to the pump must be designed to withstand the pneumatic performance of the pump

- Install the pump so that the fan can draw in sufficient cooling air
- Fit the pump at the highest point in the system, so that condensate cannot collect in the head of the pump.

SERVICE REQUIREMENTS

The diaphragm and valves are the only parts of the KNF diaphragm pumps subject to wear. They are easy to change, as no special tools are needed.

Please visit our website at <u>www.knf.com</u> for further information.

KNF NEUBERGER, INC.

2 Black Forest Road Trenton, NJ 08691-1810 609-890-8600 • Fax 609-890-8323



Diaphragm stabilization system

An additional diaphragm, the stabilization diaphragm, separates the underside of the working diaphragm from the "crankcase" of the pump (see Fig. 2). The space between the two diaphragms (called a vacuum chamber) is connected with the suction side of the pump via a balancing connection. This way, the vacuum chamber has approximately the same pressure as the working space of the diaphragm pump. The pressure difference between the upper and underside of the diaphragm approaches zero. The working diaphragm remains stable, independent of the inlet pressure of the pump. This improves the suction speed of the pump significantly, over its entire working range.



Fig. 1: Diaphragm behavior, due to the pressure difference between working space and "crank" space (normal diaphragm pump)



Fig. 2: Diaphragm stabilization system with additional diaphragm. This improves the suction speed of the pump significantly, over its entire working range.